



Are you hesitant about spinal cord experiments on rats?

Spinal Cord Clamps for Rats **STS-B**

These spinal cord clamps for rats are used in conjunction with SR series stereotaxic instruments. To ensure that the fragile spinal cord of a rat is clamped securely and correctly, this model is designed so that the clamp feeling is communicated to the operator's hand.



Installable not only on SR-5R/6R, but also on SR-5N/6N you already own



A combination of SR-5R + STS-B

Used in combination, they serve as a successor to the ST-7 which was popular for many years as a Stereotaxic Instrument with Spinal Cord Clamps for Rats.

The SR-5R comes with the SM-15 Stereotaxic Micromanipulator.

