


Double-ended Magnet System – Key Factors in the Choice of Equipment



Setup Examples

- Setup Example 1 : Manipulators are directly installed on double magnetic stands and placed beside the microscope.
- Setup Example 2 : Isolation Table System configured for brain slice patch clamp recording.
- Setup Example 3 : The stage is configured to allow manipulators to be situated in the proximity of the specimen.
- Setup Example 4 : A simple setup for brain slice patch clamp recording.
- Setup Example 5 : A stage and manipulators can be installed with any type of microscopes.
- Setup Example 6 : A setup with a stage and manipulators alone without using a microscope.